



Surgical Burs



الماس تک، همیشه تیز




Accurate high quality ALMASTECH burs specially designed to use in:

1. Neurosurgery
2. ENT surgery
3. Orthopedics
4. Dermabrasion
5. Maxilo-facial surgery

ALMASTECH burs have been classified in 5 models to meet the needs of surgeons:

- | | |
|---|---|
| 1. Fine diamond burs |  |
| 2. Coarse diamond burs |  |
| 3. Super-coarse diamond burs |  |
| 4. Ultra-coarse diamond burs |  |
| 5. Cutting burs |  |







- Highest ALMASTECH quality standards
- Instrument length of 50 – 100 mm available
- Working parts made of stainless steel and diamond
- Working part diameter from 0.6 – 7.0mm receivable
- ALMASTECH burs have been designed for working on bone, soft issue and cements

 $\emptyset = \text{mm}$	 0.8 1 2 3 4 5 6 7
Length 70 mm REF number	AT1001.1070.008 AT1001.1070.010 AT1001.1070.020 AT1001.1070.030 AT1001.1070.040 AT1001.1070.050 AT1001.1070.060 AT1001.1070.070

- Fine diamond burs are used for final polishing of bones and soft tissues
- Diamond size of this mode is about 50 micron
- Fine diamond burs are commonly used in dental, ENT and Neuro surgeries
- Stainless steel grade of these burs are 420.
- These burs are reusable and distributed in Non-Sterile state.

 $\varnothing = \text{mm}$							
Length 70 mm REF number	AT1001.1070.120	AT1001.1070.124	AT1001.1070.130	AT1001.1070.140	AT1001.1070.150	AT1001.1070.160	AT1001.1070.070
Length 100 mm REF number	---	---	AT1001.1100.130	AT1001.1100.140	AT1001.1100.150	---	---

- Coarse diamond burs are used to start the polishing procedure and also forming the soft bones
- Coarse diamond burs are commonly used in ENT, Neuro and Maxillo-facial surgeries.
- Diamond size of this model is about 180 micron
- Stainless steel grade of these bure are 420.
- These burs are reusable and distribute in Non-Sterile state.

					
Ø = mm	3	4	5	6	7
Length 70 mm REF number	AT1001.1070.230	AT1001.1070.240	AT1001.1070.250	AT1001.1070.260	AT1001.1070.270


- Super-coarse diamond burs are used to scrape the bones more accurate than cutting burs
- Coarse diamond burs are commonly used in ENT, Neuro and Maxillo-facial surgeries.
- Diamond size of this model is about 250 micron
- Stainless steel grade of these burs are 420.
- These burs are reusable and distribute in Non-Sterile state.

 $\emptyset = \text{mm}$	4.6	5.7	6.7	8
Length 70 mm REF number	AT1002.1070.346	AT1002.1070.357	AT1002.1070.367	AT1002.1070.380

- Ultra-coarse diamond burs are used to scrape the bones with different procedure besides cutting burs.
- Ultra-coarse diamond burs are commonly used in ENT, Neuro and maxillo-facial surgeries.
- Diamond size of this model is about 350 micron
- Stainless steel grade of these burs are 420.
- These burs are reusable and distribute in Non-Sterile state.



Ø = mm



0.6 1 2 2.4 3 4 5 6 7

Length 70 mm
REF number

AT1003.1070.406
AT1003.1070.410
AT1003.1070.420
AT1003.1070.424
AT1003.1070.430
AT1003.1070.440
AT1003.1070.450
AT1003.1070.460
AT1003.1070.470

- Cutting burs are commonly used in ENT, neuro and maxillo-facial surgeries.
- Stainless steel grade of these bure are 420.
- These burs are reusable and distribute in Non-Sterile state.

برخی مشتریانی که به ما اعتماد کرده اند :



 almas_tech_company

 almastech.co@gmail.com

 +98 935 452 0734

 www.almas-tech.com

 +98 910 006 8425

 Modarres Science and Technology
Park, Babaei Highway, Tehran, Iran